

## The exofin® HVTA Advantage



OVERALL LENGTH  
66mm



OVERALL LENGTH  
44mm

## exofin® HVTA Comparison to derma+flex® QS™

### EASY APPLICATION

Screw on applicator and begin applying to patient while squeezing tube to express desired amount of glue.

### HANDLING

exofin® HVTA applicator is longer and wider which makes it easier to assemble with gloved hands.

### NO FAULT APPLICATION

No special technique required. Clinician may apply adhesive as thick or thin as they desire.

### FASTER DRYING TIMES

The amount of glue applied does not change the set time of the adhesive, as long as it is applied in an even layer on the skin. Typical "dry" times are less than 30 seconds.

### VISCOSITY

exofin® high viscosity formula is the highest in the tissue adhesive category. It will not run in to unintended areas when applied correctly.

### VERY LITTLE EXOTHERM

Very little heat is generated, thus patient comfort in the ER is still a top priority for exofin® HVTA.

### EASY TRANSITION

Current derma+flex® QS™ users will see very little difference in the appearance of the packaging, but will immediately see the benefits of exofin® HVTA.

**exofin® HVTA Product Codes:** Case - EX70407-01 • Box - EX70407

To place an order, email [customerservice@chemencemedical.com](mailto:customerservice@chemencemedical.com) or call **844-MED-GLUE (844.633.4583)**

**exofin®**  
HIGH VISCOSITY TISSUE ADHESIVE

200 Technology Drive, Alpharetta, GA 30005 | **844.633.4583**  
**chemencemedical.com**

  
**CHEMENCE®**  
MEDICAL